

Cryogenic RF Switch Driver

- ✓ An easy way to drive RF switches in low-temperature cryostats
- ✓ Tuneable switching pulses
- ✓ Remote SCPI interface
- ✓ Touchscreen interface



Overview

The Cryogenic RF Switch Driver is a device to drive RF switches in cryogenic experiments. Switching pulses from a general purpose source can result in unwanted excessive heating of the experiment, resulting in a longer recovery time to reach the base temperature. The driver allows tuning the switching pulse properties to minimize the heating load. The outputs are two 8-pin connectors for easy wiring of RF switches with many ports. The driver offers a friendly local touchscreen GUI and a serial SCPI interface for remote operation.


Applications

- Cryogenic RF calibration with thruline or SOLT
- Switching between different cryogenic experiments without warm-up


Features

- Two pulse output ports, 8-pin euro-style connectors
- Each port can drive a SP6T latching RF switch
- Pulses are generated by an internal current source
- Current direction is switched to set or reset an RF port
- Pulses can be tuned in length, current, rise and fall time
- Pulse current, range: 0mA – 750mA
- Pulse length, range: 1ms – 200ms
- Rise and fall time, range: 1ms – 8ms
- Gaussian rise and fall profile
- Pulse voltage compliance: 46V
- Required power supply: 48V, 50W minimum, with earth isolation
- Pulse lines are isolated from earth
- Pulse lines are isolated from device when not in use
- Touchscreen with user-friendly GUI
- Serial SCPI interface via USB port


Application example 1: switches with common line



Application example 2: switches with actuator-series




Application example 3: two switches of different type



Application example 4: 12xSPST switches

This circuit is equivalent to application example 2, where the two SP6T switches are replaced by 6 SPST, and actuators are wired in series



SCPI interface use example with Python/Jupyter Notebook

```
In [1]: import pyvisa
rm = pyvisa.ResourceManager()
import time
```

Open

```
In [2]: rm.list_resources()
```

```
Out[2]: ('ASRL9::INSTR',)
```

```
In [3]: RFSWDriver = rm.open_resource('ASRL9::INSTR')
RFSWDriver.write_termination = '\n'
RFSWDriver.read_termination = '\n'
RFSWDriver.baud_rate = 9600
```

```
In [4]: RFSWDriver.query('*IDN?')
```

```
Out[4]: 'RF Switch driver\r'
```

Set / reset

```
In [5]: # Set switch. Syntax: Switch (1 or 2), Port (1 to 6)
RFSWDriver.query('SET 1,1')
```

```
Out[5]: 'Sw1, P1 set.\r'
```

```
In [6]: # Reset switch. Syntax: Switch (1 or 2), Port (1 to 6)
RFSWDriver.query('RST 1,1')
```

```
Out[6]: 'Sw1, P1 reset.\r'
```

```
In [7]: # Reset ALL. Syntax: Switch (1 or 2)
RFSWDriver.query('RSTALL 2')
```

```
Out[7]: 'Sw2 reset all\r'
```

Switch status

```
In [16]: # Port Status
RFSWDriver.query('STATUS')
```

```
Out[16]: 'SW1: P1:R || P2:S || P3:S || P4:S || P5:R || P6:R ||| SW2: P1:R || P2:R || P3:R || P4:S || P5:S || P6:S \r'
```

Read or set pulse parameters

```
In [17]: # Set Rise Time. Syntax: Switch (1 or 2), Port (1 to 7) or 0, RiseTime ms (0 to 8).
#If Port=0, all the ports are set the same
#If Port=7, the setting is for Reset ALL
RFSWDriver.query('RISETIME 1,1, 2')
```

```
Out[17]: 'Rise Time updated\r'
```

```
In [18]: # Read Rise Time. Syntax: Switch (1 or 2), Port (1 to 7)
RFSWDriver.query('RISETIME 1,1')
```

```
Out[18]: '2\r'
```

```
In [19]: # Set Fall Time. Syntax: Switch (1 or 2), Port (1 to 7), Fall Time ms (0 to 8)
RFSWDriver.query('FALLTIME 1,1,5')
```

```
Out[19]: 'Fall Time updated\r'
```

```
In [20]: # Read Fall Time. Syntax: Switch (1 or 2), Port (1 to 7)
RFSWDriver.query('FALLTIME 1,1')
```

```
Out[20]: '5\r'
```

```
In [21]: # Set Length. Syntax: Switch (1 or 2), Port (1 to 7), pulse length ms (0 to 1000)
RFSWDriver.query('LENGTH 1,1,20')
```

```
Out[21]: 'Pulse length updated\r'
```

```
In [22]: # Read Length. Syntax: Switch (1 or 2), Port (1 to 7), pulse length ms (0 to 1000)
RFSWDriver.query('LENGTH 1,1')
```

```
Out[22]: '20\r'
```

```
In [23]: # Set pulse current. Syntax: Switch (1 or 2), Port (1 to 7), pulse current mA (0 to 750)
RFSWDriver.query('CURR 1,1,400')
```

```
Out[23]: 'Pulse current updated\r'
```

```
In [24]: # read pulse current. Syntax: Switch (1 or 2), Port (1 to 7)
RFSWDriver.query('CURR 1,1')
```

```
Out[24]: '400\r'
```

```
In [25]: # Preset
RFSWDriver.query('PRESET')
```

```
Out[25]: 'Preset\r'
```

Pictures



Power and USB ports




Pulse ports, 2x8-pin euro-style



Main menu, with toggling buttons



Standby mode



Pulse library

Notes

- The driver can work only with latching relays
- The driver does not actively check the RF port status or the success of a switching event
- The port status displayed corresponds to the last pulse sent
- Pulse lines are isolated from earth only if an isolated power supply is used